

# Chemicals used in the lab provoke nauseas, headache, tiredness, vertigo


The health of your chemists  
is in **DANGER!**



**\* NEW!**

Upgrade **immediately** their  
protection with filtration  
fume hoods and storage  
cabinets **captair®**

- Excellent filtration efficiency according to NFX 15-211 standard
- No duct to outside, a simple electrical socket is sufficient to connect a captair®
- Mobile, can be moved anytime to the most suitable location
- Energy saving: running costs 8-10 times cheaper than traditional Ducted Fume Hood



**captair® flex™**  
by erlab

Ductless Filtering Fume Hood

\* equipped with  
new flexible filtration columns to  
filter solvents, acids, bases, powders.  
Also suitable for use in clean rooms

**captair® store**  
by erlab

Ductless Filtering Storage Cabinet

[www.captair.com](http://www.captair.com)

Please return by Fax to: +604 397 8221

DR. RAHMATULLAH HOLDINGS SDN. BHD.  
(478017-T)

6716 Pengkalan Tambang,  
Permatang Pasir, 13500 Bukit Mertajam,  
Pulau Pinang.


Email: mamagani1@yahoo.com  
Tel: 04-3901971 Fax: 04-397 8221

You are exposed to a significant health risk if you inhale every day more than the Threshold Limit Value (TLV) of any chemical.

For example, TLVs of some chemicals:

Toluene:	50 ppm
Trichlo:	75 ppm
Nitric acid:	2 ppm
Sodium hydroxide	2 mg/m <sup>3</sup>

The TLVs valid in the USA, called TWA (Time Weight Average), are determined by the ACGIH (American Conference of Governmental Industrial Hygienists). We have the pleasure to offer you a free Poster listing the currently available TWAs. Just fill out and return the coupon on the back to us!



A green product from **erlab®**  
group

# captair® safety enclosures

FOR THE SAFE HANDLING OF YOUR HAZARDOUS AND SENSITIVE PRODUCTS

## captair® bio by erlab

biocap™ PCR workstations to protect the DNA sample against DNAs/RNases contamination from the operator, from the room and from the previous experiment!

## captair® flow by erlab

flowcap™ vertical laminar flow cabinets are designed to perform handling of non-pathogenic materials in an ultra-clean environment: cell culture, in vitro culture, microbiology, electronics, optic, etc. . .

## captair® pyramid by erlab

pyramid™ 2200 disposable glove box for the protection of the user and sensitive products, i. e. opening of suspicious packages, protection of oxidizable products, fingerprint analysis, etc. . .

## captair® flex™ by erlab

Ductless filtering fume hoods with modular filtration columns

For the handling and weighing of liquid chemicals and powders, totally eliminated by filtration.

A range of 17 models, all in compliance with NF X15-211 standard. Mobile, energy saving and environmental friendly!

## captair® store by erlab

Ductless filtering storage cabinets

Working 24 hours a day, captair® store cabinets help not only to clean the residual vapors released from the bottles stored inside, but they will also clean permanently the air of the laboratory room, providing full protection from chemicals, toxics and odours.

## shelf™, ministore™, AVP™

models: For the storage of chemicals, totally eliminated by filtration. A range of 11 models!

[www.captair.com](http://www.captair.com)

I would like to receive information on:

captair® flex  
Ductless filtering fume hoods

captair® store  
Ductless filtering storage cabinets

captair® bio  
biocap™ PCR workstations

captair® flow  
Vertical laminar flow workstations

captair® pyramid®  
Disposable glove box

My applications involve:

- Handling  Weighing  Storage of:  
 Acids  Solvents  Bases  Powders in:  
 Standard chemistry laboratory  Clean room laboratory class.....

Others: \_\_\_\_\_

Please send me your free TLV poster

Please send me a costs comparison sheet to help me comparing the energy costs between a Ducted Fume Hood and a captair® filtering fume hood.



Please complete your information:

Your name: \_\_\_\_\_

Your position: \_\_\_\_\_

Organisation name: \_\_\_\_\_

Address: \_\_\_\_\_

\_\_\_\_\_

Tel. : \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

Please return the completed coupon by fax to: +604 397 8221

Alternatively, return it by mail to: **mamagani1@yahoo.com**



**DR. RAHMATULLAH HOLDINGS SDN. BHD.**  
 6716 Pengkalan Tambang  
 Permatang Pasir  
 13500 Bukit Mertajam  
 Phone: +604 390 1971